


BALTON[®]

www.balton.pl



Dialysis

DIALYSIS

DIALYSIS SET

| | |
|----------------------------------|---|
| Single lumen, straight | 2 |
| Single lumen "Y" | 2 |
| Double lumen | 3 |
| Triple lumen | 4 |

DIALYSIS CATHETER

| | |
|----------------------------|---|
| Single lumen | 5 |
| Single lumen "Y" | 5 |
| Double lumen | 6 |
| Triple lumen | 6 |

HEMODIALYSIS NEEDLE & ACCESSORIES

| | |
|---|---|
| Hemodialysis needle | 7 |
| Dialysis extender for dialyzer rinsing / cleaning | 8 |
| Dialysis extender "Y" | 8 |

DIALYSIS SET single lumen, straight



Set components:

single lumen straight catheter
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

| Reference: guide wire "J" type | Reference: guide wire Straight type | Catheter size / length | Dilator I / II | Guide wire size / length |
|--------------------------------------|---|---------------------------|-------------------|-----------------------------|
| ZDJP8F15 | ZDJP8F15P | 8F / 15 cm | 8F /10F | .035" / 60 cm |
| ZDJP8F18 | ZDJP8F18P | 8F / 18 cm | 8F /10F | .035" / 60 cm |
| ZDJP8F20 | ZDJP8F20P | 8F / 20 cm | 8F /10F | .035" / 60 cm |

When ordering use Ref. Number, e.g.: ZDJP8F20.

NOTE: On request we provide catheters with hydrophilic surface protection. When ordering use Ref. Number, e.g.: ZDJP8F20H.

DIALYSIS SET single lumen "Y"



Set components:

single lumen catheter "Y"
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

| Reference: guide wire "J" type | Reference: guide wire Straight type | Catheter size / length | Dilator I / II | Guide wire size / length |
|--------------------------------------|---|---------------------------|-------------------|-----------------------------|
| ZDJY8F15 | ZDJY8F15P | 8F / 15 cm | 8F /10F | .035" / 60 cm |
| ZDJY8F18 | ZDJY8F18P | 8F / 18 cm | 8F /10F | .035" / 60 cm |
| ZDJY8F20 | ZDJY8F20P | 8F / 20 cm | 8F /10F | .035" / 60 cm |


When ordering use Ref. Number, e.g.: ZDJY8F20.

NOTE: On request we provide catheters with hydrophilic surface protection. When ordering use Ref. Number, e.g.: ZDJY8F20H.

DIALYSIS SET double lumen

Set components:

double lumen catheter (Straight type or Curved type)
 needle 18G, length 7 cm
 for catheter size 6.5F - needle 19G, length 5 cm
 guide wire
 dilatator I
 dilatator II



| Reference: guide wire "J" type | Reference: guide wire Straight type | Catheter size / length | Dilatator I / II | Guide wire size / length | Version |
|--------------------------------------|---|---------------------------|---------------------|-----------------------------|-----------|
| ZDD6.5F10 | ZDD6.5F10P | 6.5F / 10 cm | 7F / 8F | .027" / 60 cm | PEDIATRIC |
| ZDD6.5F12 | ZDD6.5F12P | 6.5F / 12 cm | 7F / 8F | .027" / 60 cm | |
| ZDD8F10 | ZDD8F10P | 8F / 10 cm | 8F / 10F | .035" / 60 cm | |
| ZDD8F12 | ZDD8F12P | 8F / 12 cm | 8F / 10F | .035" / 60 cm | |
| ZDD8F15 | ZDD8F15P | 8F / 15 cm | 8F / 10F | .035" / 60 cm | |
| ZDD11F15 | ZDD11F15P | 11F / 15 cm | 10F / 12F | .035" / 70 cm | |
| ZDD11F18 | ZDD11F18P | 11F / 18 cm | 10F / 12F | .035" / 70 cm | |
| ZDD11F20 | ZDD11F20P | 11F / 20 cm | 10F / 12F | .035" / 70 cm | |
| ZDD12F15 | ZDD12F15P | 12F / 15 cm | 10F / 12F | .035" / 70 cm | |
| ZDD12F18 | ZDD12F18P | 12F / 18 cm | 10F / 12F | .035" / 70 cm | |
| ZDD12F20 | ZDD12F20P | 12F / 20 cm | 10F / 12F | .035" / 70 cm | |

When ordering use Ref. Number, e.g.: ZDD12F20.

NOTE: other configurations available


H - catheter with hydrophilic coating

When ordering use Ref. Number, e.g.: ZDD11F20HP.

Z - catheter with curved arms

When ordering use Ref. Number, e.g.: ZDD11F20ZP.

DIALYSIS SET triple lumen



Set components:

triple lumen catheter (Straight type or Curved type)
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

| Reference: | Catheter size / length | Dilatator I / II | Guide wire size / length |
|------------|------------------------|------------------|--------------------------|
| ZDT12F15 | 12F / 15 cm | 10F / 12F | .038" / 70 cm |
| ZDT12F18 | 12F / 18 cm | 10F / 12F | .038" / 70 cm |
| ZDT12F20 | 12F / 20 cm | 10F / 12F | .038" / 70 cm |

When ordering use Ref. Number, e.g.: ZDT12F15.

NOTE: other configurations available:

H - catheter with hydrophilic coating

When ordering use Ref. Number, e.g.: ZDT12F15H

Z - catheter with curved arms

When ordering use Ref. Number, e.g.: ZDT12F15Z.

DIALYSIS CATHETER single lumen



Reference:

KDJP8F15
KDJP8F18
KDJP8F20

Catheter size / length

8F / 15 cm
8F / 18 cm
8F / 20 cm

When ordering use Ref. Number, e.g.: KDJP8F20.

NOTE: Catheter with curved arms available to order. When ordering use Ref. Number, e.g.: KDJP8F15Z.

DIALYSIS CATHETER single lumen "Y"



Reference:

KDJY8F15
KDJY8F18
KDJY8F20

Catheter size / length

8F / 15 cm
8F / 18 cm
8F / 20 cm

When ordering use Ref. Number, e.g.: KDJY8F20.

DIALYSIS CATHETER double lumen



| Reference: | Catheter size / length | Version |
|------------|------------------------|------------|
| KDD6.5F10 | 6.5F / 10 cm | PEDIATRIC |
| KDD6.5F12 | 6.5F / 12 cm | |
| KDD8F10 | 8F / 10 cm | |
| KDD8F12 | 8F / 12 cm | |
| KDD8F15 | 8F / 15 cm | |
| KDD11F15 | 11F / 15 cm | FOR ADULTS |
| KDD11F18 | 11F / 18 cm | |
| KDD11F20 | 11F / 20 cm | |
| KDD12F15 | 12F / 15 cm | |
| KDD12F18 | 12F / 18 cm | |
| KDD12F20 | 12F / 20 cm | |

When ordering use Ref. Number, e.g.: KDD8F12.

NOTE: Catheter with curved arms available to order. When ordering use Ref. Number, e.g.: KDD8F15Z.
H - catheter with hydrophilic coating. When ordering use Ref. Number, e.g.: KDD12F15H.

DIALYSIS CATHETER triple lumen



| Reference: | Catheter size / length |
|------------|------------------------|
| KDT12F15 | 12F / 15 cm |
| KDT12F18 | 12F / 18 cm |
| KDT12F20 | 12F / 20 cm |

When ordering use Ref. Number, e.g.: KDT12F15.

NOTE: other configurations available
H - catheter with hydrophilic coating. When ordering use Ref. Number, e.g.: KDT12F15H
Z - catheter with curved arms. When ordering use Ref. Number, e.g.: KDT12F15Z.

HEMODIALYSIS NEEDLE

The needle is made of special steel alloy that ensures exceptional smoothness of its surface. Thin walls and a triple - grinded spike punctures the fistula delicately and precisely. The surface of the needle is covered with silicone and therefore it goes through the tissues easily and well tolerated in the vessel. Specially shaped wings hold the needle easily and allow easy hanging. Moving wings make it easy to attach the needle to the skin - the position of the butterfly - nut can be changed without turning the needle in the vessel. The clamps that are easy to fasten and reliable in action are colored according to the line concerned: blue in the case of venous needles, red in the case of arterial needles, that allows easy application and prevents mistakes.



| Reference: | Needle size | Version | Wings colour | Clamps colour |
|------------|-------------|----------|--------------|---------------|
| ITD1.5LL | 17G | Arterial | red | red |
| ITD1.6LL | 16G | Arterial | green | red |
| ITD1.8LL | 15G | Arterial | yellow | red |
| IZD1.5LL | 17G | Venous | red | blue |
| IZD1.6LL | 16G | Venous | green | blue |
| IZD1.8LL | 15G | Venous | yellow | blue |
| IDY1.5LL | 17G | "Y" | red | red/blue |
| IDY1.6LL | 16G | "Y" | green | red/blue |
| IDY1.8LL | 15G | "Y" | yellow | red/blue |

When ordering use Ref. Number, e.g.: IZD1.5LL.

DIALYSIS EXTENDER for dializator rinsing / cleaning


Reference:

PD7

Length

7 cm

Diameter

4.6 / 3.3 mm

When ordering use Ref. Number, e.g.: PD7.

DIALYSIS EXTENDER "Y"




"Y" extender is used to transform straight catheter into a "Y" catheter.

Reference:

PUDY

When ordering use Ref. Number, e.g.: PUDY.



Producer of Medical Equipment

294 Modlińska Str., 03-152 Warsaw, POLAND
ph.: (+48) 22 597 44 00, fax: (+48) 22 597 44 44
e-mail: balton@balton.pl, www.balton.pl